GlutaMAX[™]-I (100X)

Description

GlutaMAX[™]-I, L-alanyl-L-glutamine, is a dipeptide substitute for L-glutamine. GlutaMAX[™]-I can be used as a direct substitute for L-glutamine at equimolar concentrations in both adherent and suspension mammalian cell cultures with minimal or no adaptation. GlutaMAX[™]-I is highly soluble, heat-stable, and improves growth efficiency and performance of mammalian cell culture systems. GlutaMAX[™]-I eliminates problems associated with the spontaneous breakdown of L-glutamine into ammonia during incubation, allowing for longer lasting cultures. GlutaMAX[™]-I supplement is also available pre-supplemented in a variety of media formulations.

Product	Catalog No.	Amount	Storage	Shelf Life*	
GlutaMAX [™] -I	35050-061 35050-079 35050-038** 35050-087**	100 mL 20 × 100 mL 100 mL 20 × 100 mL	15°C to 30°C	24 months	

* Shelf Life duration is determined from Date of Manufacture.

** For European customers only.

Intended Use

For in vitro diagnostic use.

Important Information

• Formulated as a 200 mM (100X) L-alanyl-L-glutamine dipeptide in 0.85% NaCl

Safety Information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Use

Direct equimolar substitution for L-glutamine in mammalian cell culture.

Supplement with GlutaMAX[™]-I

- 1. Aseptically add the required molar concentration to culture medium.
- If a lag in growth is observed upon conversion to GlutaMAX[™]-I, it is recommended to supplement the culture media with GlutaMAX[™]-I and L-glutamine at a 3:1 molar ratio (e.g., for 8 mM final concentration, use 6 mM GlutaMAX[™]-I and 2 mM L-glutamine) for the first and second passages during adaptation to the new culture conditions.

Related Products

Product	Catalog No.	
L-glutamine 200mM (100X), liquid	25030	
GlutaMAX [™] -I CTS [™]	A12860	
DMEM/F-12, GlutaMAX [™]	10565	
RPMI, GlutaMAX™	61870	
Opti-MEM [®] Reduced Serum Medium, GlutaMAX [™]	51985	

Explanation of Symbols and Warnings

The symbols present on the product label are explained below:

<u>/</u> !	\triangle		IVD		STERILE A		i	
Caution, accompanying		In vitro diagnostic medical device		Sterilized using aseptic processing techniques				
	REF	***		LOT	Ce	-	X	
Use By:	Catalog number	Manufacturer	Ba	atch Code	Europeai Communi		Temperature Limitation	

Limited Use Label License

No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact **outlicensing@lifetech.com** or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **www.lifetechnologies.com/termsandconditions**. If you have any questions, please contact Life Technologies at **www.lifetechnologies.com/support**.

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit www.lifetechnologies.com/support For further assistance, email **techsupport@lifetech.com**

© 2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

